

**SOP:** Propagation of human renal cortical epithelial cells (HRCEpiC, Sciencell)  
**Date modified:** 6/03/09  
**Modified by:** A. Haydock/M. Dorschner (UW)

### **Ordering Information**

Human renal cortical epithelial cells (HRCEpiC) may be ordered as frozen ampoules. Each vial contains  $>5 \times 10^5$  cells in a 1ml volume.

To order frozen ampoules + media:

Name:	HRCEpiC – Human Renal Cortical Epithelial Cells
Item #:	4110 (HRCEpiC - Cryopreserved ampoule)
	4101 (EpiCM, Epithelial Cell Medium)

### **Notes:**

The number of media kits purchased depends on the target number of cells to be generated. A rule of thumb is 10 media kits for every initial cryopreserved ampoule. It is strongly recommended to purchase all of the media that will be required for a complete expansion series, since media supply may be erratic. All culture flasks must be coated with poly-L-lysine ( $2\mu\text{g}/\text{cm}^2$ ).

### **Materials List**

1. Cell-type specific medium (Sciencell)
2. Poly-L-Lysine (#0413, 10mg/ml, Sciencell)
3. T225 culture flasks
4. Graduated pipets (1, 5, 25mL)
5. Pen-strep solution (if required; Lonza typically supplies antibiotics)
6. Hemocytometer
7. Micropipet w/ P20 tips
8. Microscope

### **Procedure**

#### **A. Initiation of culture from cryopreserved cells**

- 1) Rapidly thaw cells by holding vial and gently rotating in  $37^\circ\text{C}$  water bath.
- 2) As soon as ice crystals disappear, dispense contents of vial into a flask at  $7500 \text{ cells}/\text{cm}^2$  density.
- 3) Let cells recover for 16 hours in a  $37^\circ\text{C}$ , 5%  $\text{CO}_2$  humidified incubator.

#### **B. Sub-culture**

- 1) Propagate cells until density reaches 70-80% confluence.
- 2) Decant medium.
- 3) Wash cells with warm 1X PBS.
- 4) Add 10mLs of Trypsin and return to incubator for 10-15 minutes.
- 5) Immediately remove cells and pellet at  $500 \text{ xg}$  for 5 minutes ( $4^\circ\text{C}$ )
- 6) Wash cells 2X with 1X PBS.

- 7) Gently re-suspend cell pellet in warm medium.
- 8) Count cells with hemocytometer.
- 9) Add warmed medium to flasks.
- 10) Seed flasks at **7,500 cells/cm<sup>2</sup>**
- 11) Record each subculture event as a passage

### **C. Maintenance**

- 1) Change media the day after seeding and every OTHER day thereafter.
- 2) Increase media volume as confluency increases (volumes assume the use of
- 3) T225 flasks):
  - a. 25 % = 1mL/5 cm<sup>2</sup>
  - b. 25-45% = 1.5mL/ 5 cm<sup>2</sup>
  - c. 45%+ = 2mL/ 5 cm<sup>2</sup>.
- 4) Per the above an exemplary schedule might be:
  - a. day 1, plate into T225: use 50 mls of media.
  - b. day 2, change media, use 50 mls of media
  - c. day 4, change media, use 100 mls of media (if confluency is >50%)
  - d. day 6, change media, use 100 mls of media (or harvest if ready)
  - e. day 7 or 8 (harvest when cells reach  $6 \times 10^6$  cells/flask)

### **D. Harvest**

- 1) Pass cells 3-4 times until the desired cell number is achieved (primary cells will senesce after 4-5 passages).
- 2) Remove cells from flasks according to protocol described above under 'subculturing'
- 3) Examine viability using trypan blue staining (SOP TP-7)